

Roho Comparison Chart

| Select Air System Specifications: | <i>SelectAir</i> | <i>SelectAir Max</i> | <i>BariSelect</i> (all sizes) | <i>Hybrid Select</i> |
|--|----------------------|----------------------|----------------------------------|----------------------|
| Mattress System Type | Mattress replacement | Mattress replacement | Mattress replacement | Mattress Overlay |
| Patient Weight Limit | 400 lbs. | 400 lbs. | 680 lbs. | 300 lbs. |
| One Year Limited Warranty | X | X | X | X * |
| Control Unit Specifications: | | | | |
| UL 544 Certification # E183348 | X | X | X | X |
| Current - Less than 3 amps | x | x | x | x |
| Voltage - 120V AC 50/60 Hz | x | x | x | x |
| Flow Rate 80-100 SLPM | x | x | x | x |
| Mattress, Cover, Foam, Sheet Specifications | | | | |
| Flammability - 16 CFR Part 1632 (composite testing) | x | x | x | x |
| Latex free | x | x | x | x |
| Surface Resistivity - Anti Static | x | x | x | x |
| Bacterial Resistance - Meets AATCC-147-1993 | x | x | x | x |
| Control Unit Features: | <i>SelectAir</i> | <i>SelectAir Max</i> | <i>BariSelect</i> (all sizes) | <i>Hybrid Select</i> |
| Softer/Firmer - Push to increase (firmer) or decrease (softer) comfort level. | x | x | x | x |
| Upright - Increases pressure to compensate for the upright patient. | x | x | x | x |
| Firm - Fully inflates mattress to aid in placement and maneuvering of patient. | x | x | x | x |
| System Lock - Locks the system to deter unauthorized personnel from altering the settings. | | x | x | x |
| Auto Comfort - Automatically adjusts the system according to the patient's weight distribution. | | x | x | |

Roho Comparison Chart

| Control Unit Features (con't): | SelectAir | SelectAir Max | BariSelect (all sizes) | Hybrid Select |
|---|-----------|---------------|------------------------|---------------|
| Pulsation - Continually readjusts system pressure. SelectAir Max & BariSelect -Low-every 15 minutes, Medium-every 7 minutes, High-every minute. HybridSelect - Low-every 5 minutes, High-every minute. | | x | x | x |
| CPR/Emergency Deflate - When engaged, allows quick deflation of the system in case of an emergency. | x | x | x | x |
| Digital Display - (LCD) - Easy to read digital display. | | | x | |
| Remote Capability - An optional hand control is available (additional cost) which allows the patient to self-adjust the softer/firmer and power features. | | | x | |
| Back-lit Keys - Allow system features to be visible in dim/dark environments. | | | x | |